
## Standard Temperature Data Logger -90 $^\circ\text{C}$ to $\,+\,140$ $^\circ\text{C}$



- + Collects and stores temperature information in the internal memory
- + Recover the data through the docking station using XpertLog Software
  - ITS-90 calibration coefficients are stored in the internal memory
- + Complies with FDA 21 CFR Part 11

+

+ Network ready: multiple users, centralized, server database, active directory authentication





Visit www.lives-international.com for more details.

----->

XpertLog software ----->



LVSC.3.1.8 - May 2020